


This is how much power each electrode sends

What are intensities?

into your muscles. In the app, these appear as percentages (0-100%) next to each muscle group, and are controlled by the + / - buttons.

10:41



If you lift tiny weights every day, you'll stay active but you won't get stronger. EMS works the same way: if your intensities are too low, you won't see real results.

Here's how to optimize your intensities:

Train at Your Highest Tolerable Intensity

right at the edge where you can still complete movements

This is the #1 rule.

with proper form - but just barely! That's your Highest Tolerable Intensity. You should feel your muscles working hard, not just tingling. Increase Intensity Within Each Session

Push yourself to where the impulses feel really challenging -

so continuously adjust intensity upward. What feels challenging at 15% intensity in minute 2 will feel too easy by minute 5. Keep

increasing throughout the entire 20-minute session.

EMS becomes easier as your muscles adapt during the session,


Rule of Thumb: Aim to roughly double your starting intensity by session's end. Example:

Larger muscles tolerate much higher stimulation than smaller

Start: 15% intensity (minute 2)

Finish: 30% intensity (minute 19)

Pro Tip:



ones. It's normal to reach 50%

intensity on legs and glutes while

keeping arms and chest around

20% in the same session.

Increase Intensity Over Time Like progressive overload at the gym, you need to keep challenging your muscles.

4 Use Soreness as Feedback

intensities that are too conservative.

Remember:

Muscle soreness is normal after EMS, especially at first. If you

feel no soreness or fatigue, you're probably training at

Intensity Beats Frequency Two high-effort sessions per week will

deliver better results than four easy ones. Push yourself, track your progress, and keep dialing it up — your results will follow.



0

Terms and Conditions. You can unsubscribe at any time.

marketing communication at sqai.co. For more information, please refer to our Privacy Policy and



11900 Biscayne Blvd Ste 105 Miami, FL 33181